



SCAN▶X® Duo

SCAN▶X® Duo

Intraoral Digital Imaging

AIR
TECHNIQUES equipped for life®

SEE MORE.

ScanX Duo captures high resolution, detail-rich images for improved diagnostics.

- Scan up to 22 line pairs per mm and see unsurpassed digital clarity.
- Enlarge, invert, or sharpen your radiographs to maximum diagnostic information and enrich patient discussions about treatment plans.
- Automatically store and file electronic images for future use or transmittal.

DO MORE.

ScanX Duo streamlines procedures and increases office productivity.

- Work more efficiently with no workflow disruptions.
- Utilizing In-Line Erase, Phosphor Storage Plates (PSPs) are automatically scanned and erased in one continuous cycle- intraoral images and viewable in seconds.
- Improve patient comfort with thin, flexible, wireless, reusable phosphor plates in intraoral sizes 0, 1, 2, 3 and 4

SPEND LESS.

ScanX Duo, both practical and cost-effective, will pay for itself in less than a year.

- Eliminate the darkroom, film, hazardous chemicals, and expensive ongoing processor maintenance costs.
- Works seamlessly with your existing X-ray equipment.
- Reuse inexpensive Phosphor Storage Plates hundreds of times.

SPECIFICATIONS:
Part Number: D1000F
Dimensions (DxWxD): 10.375" x 10" x 12" (26.352cm x 25.4cm x 30.48cm)
Weight: 21 lbs. (9.5 Kg)
Electrical: 100 - 240V AC, 50/60 Hz
Resolution: 20 LP/mm (1000 dpi). Depending on scanning model/pixel size selected.
Theoretical Resolution: 40 LP/mm (2000 dpi)
Ports: USB
Laser Classification: Class I laser Product, Compliance with FDA 21 CFR 1040.10 and IEC 60825-1
Imaging Medium: Two processing slots accept Phosphor Storage Plates. Intraoral sizes 0, 1, 2, 3 and 4.
Accessories Included: Power Cord, Computer Connector (USB), 2 each size 2 PSP Guides, 1 each size 0 and 1 PSP Guides, 8 each size 2 PSPs; 1 Box Barrier Envelope Size 2 (Qty. 300); 1 Sample ScanX Cleaning Sheet; 1 Sample PSP Cleaning Wipe; CD with Drivers, Utilities and Instillation Manual.
Warranty: ScanX Duo is warranted to be free from defects for a period of 24 months through an authorized Air Techniques dealer.

COMPUTER RECOMMENDATIONS:
CPU SPEED: Intel Pentium IV compatible, 1.4 GHz or faster required. 3.0 GHz Intel CORE2 recommended.
Operating System(s): Microsoft Windows 7 Professional, Enterprise, or Ultimate with Service Pack 1 or later for an Intel 32-bit or an Intel 64-bit extended (x64) processor; Microsoft Windows 8.1 Professional or Enterprise for an Intel 64-bit extended (x64) processor; or Microsoft Windows 10 Professional or Enterprise for an Intel 64-bit extended (x64) processor.
RAM: 2 GB
Hard Drive: 120 GB (2GB free for image acquisition, the rest for program files and image storage.)
Video Display Adapter: 256 MB RAM
USB Port: 2.0

SCAN▶X® DUO



CORPORATE HEADQUARTERS
1295 Walt Whitman Road
Melville, New York 11747-3062
800-247-8324



ScanX is a registered trademark of Air Techniques, Inc.
© Copyright 2018
PN 9922-333 RevF

